

Examples for Control Sera with some information

Product Name	Analytes			Details	
Tumor Marker Control	AFP CA 15-3 CA 72-4 PSA TPS	CA 19-9 CEA Ferritin SCC	CA 125 Cyfra 21-1 HCG TPA	Form: Matrix:	liquid Human Serum
				Shelf Life: Stability: Shipment:	5 years at -20°C 5 years at 4°C ambient Temp.
Allergy Control	d 01 g 06 f 13	d 02 g 12 f 17	e 01 f 02 w 06	Form: Matrix:	liquid Human Serum
				Shelf Life: Stability: Shipment:	5 years at -20°C 6 months at 4°C ambient Temp.
Plasma Protein Control	Albumin Apo A1* C 3 CRP Haptoglobin IgE α2-Macroglob. κ-Light Chains Rheumat. Factor	α1-Antitrypsin Apo B C 4 Ferritin Hemopexin IgG β2-Microglob. Prealbumin Transferrin	Antithrombin III ASLO Coeruloplasmin* α1-Glycoprotein IgA IgM λ-Light Chains Retinol Bdg Prot. Transferrin Rec.	Form: Matrix:	liquid Human Serum
				Shelf Life: Stability: Shipment:	5 years at -20°C 9 months at 4°C (*1-3 months) ambient Temp.
CSF (Cerebrospinal Fluid) Control	Albumin IgA Lactate Sodium	Chloride IgG Potassium	Glucose IgM Protein	Form: Matrix:	liquid Human Serum
				Shelf Life: Stability: Shipment:	5 years at -20°C 12 months at 4°C ambient Temp.
Cardiac Marker Control	CK-MB NT-proBNP	Homocysteine Troponin I	Myoglobin Troponin T	Form: Matrix:	liquid Human Serum
				Shelf Life: Stability: Shipment:	5 years at -20°C 4-30 days at 4°C Dry Ice
Thyroid Hormone Control	Thyroglobulin free T3	TSH total T4	total T3 free T4	Form: Matrix:	liquid Human Serum
				Shelf Life: Stability: Shipment:	5 years at -20°C 6 months at 4°C ambient Temp.
Hormone Control	Aldosterone Cortisol Estriol HCG Hydroxyprogest. LH Testosterone free T4	AFP DHEA-Sulfate Ferritin β-HCG IgE Progesterone TSH Vitamin B12	CK-MB Estradiol FSH hGH Insulin Prolactin free T3	Form: Matrix:	liquid Human Serum
				Shelf Life: Stability: Shipment:	5 years at -20°C 6 months at 4°C ambient Temp.

trust and verify



BIOREF[®]

Dammweg 1 | 63776 Moembris/Germany | Tel. +49 6029 -5121 | Fax -5594 | info@bioref.com

www.bioref.com