

# Examples for Control Sera with some information

Product Name	Analytes			Details
<b>Tumor Marker Control</b> 4 x 3 ml	AFP CA 15-3 CA 72-4 Ferritin SCC	CASA CA 19-9 CEA HCG TPA	CA 125 CA 50 Cyfra 21-1 PSA TPS	Form: liquid Matrix: Human Serum Levels: 2 Shelf Life: 5 years at 4°C Stability: 5 years at 4°C Shipment: ambient Temp.
<b>Allergy Control</b> 4 x 2 ml	d 01 g 06 f 13	d 02 g 12 f 17	e 01 f 02 w 06	Form: liquid Matrix: Human Serum Levels: 2 Shelf Life: 5 years at -20°C Stability: 6 months at 4°C Shipment: ambient Temp.
<b>Plasma Protein Control</b> 3 x 2 ml	Albumin Apo A1* C 3 CRP Haptoglobin IgE α2-Macroglob. κ-Light Chains Rheumat. Factor	α1-Antitrypsin Apo B C 4 Ferritin Hemopexin IgG β2-Microglob. Prealbumin Transferrin	Antithrombin III ASLO Coeruloplasmin* α1-Glycoprotein IgA IgM λ-Light Chains Retinol Bdg Prot. Transferrin Rec.	Form: liquid Matrix: Human Serum Levels: 3 Shelf Life: 5 years at -20°C Stability: 9 months at 4°C (*1-3 months) Shipment: ambient Temp.
<b>CSF (Cerebrospinal Fluid) Control</b> 3 x 2 ml	Albumin IgA Lactate Sodium	Chloride IgG Potassium	Glucose IgM Protein	Form: liquid Matrix: Human Serum Levels: 3 Shelf Life: 5 years at -20°C Stability: 12 months at 4°C Shipment: ambient Temp.
<b>Cardiac Marker Control</b> 4 x 2 ml	CK-MB NT-proBNP	Homocysteine Troponin I	Myoglobin Troponin T	Form: liquid Matrix: Human Serum Levels: 2 Shelf Life: 5 years at -20°C Stability: 4-30 days at 4°C Shipment: Dry Ice
<b>Thyroid Hormone Control</b> 3 x 3 ml	Thyroglobulin free T3	TSH total T4	total T3 free T4	Form: liquid Matrix: Human Serum Levels: 3 Shelf Life: 5 years at -20°C Stability: 6 months at 4°C Shipment: ambient Temp.
<b>Hormone Control</b> 3 x 6 ml	Aldosterone Cortisol Estrinol HCG Hydroxyprogest. LH Testosterone free T3 Vitamin B12	AFP DHEA-Sulfate Ferritin β-HCG IgE Progesterone TSH total T4	CK-MB Estradiol FSH hGH Insulin Prolactin total T3 free T4	Form: liquid Matrix: Human Serum Levels: 3 Shelf Life: 5 years at -20°C Stability: 6 months at 4°C Shipment: ambient Temp.

trust and verify



**BIOREF**<sup>®</sup>

Dammweg 1 | 63776 Moembris/Germany | Tel. +49 6029 -5121 | Fax -5594 | info@bioref.com

www.bioref.com